

MK-A07

Accompany Chair



*Better care
Better comfort!*

MK-A07

ACCOMPANY CHAIR

DESCRIPTION

MK-A07 accompany chair is for the patient's attendant in hospital, it can be fold as a chair in day time, while unfold as a bed in night, which saves space and is safe and reliable. The elegant arc-shaped armrests are ergonomically designed, comfortable and beautiful, and flexible to use.

FEATURES

- Accompany chair can be open as a simple hospital bed, it can be folded into chair.
- The structure is made of powder coating steel.
- Frame available, PVC and high quality sponge applied for chair surface.
- The accompany chair is equipped with Inside extension spring.
- Suitable for high-class patient room or office.
- Two pieces 3 inch noiseless wheels, four small wheels.

STANDARD SPECIFICATION

Folding Size	740*620*350mm
Unfolding Size	1900*620*580mm
Gross Weight	31KG
Packing Size	740*660*380mm



Seat extension

seat can be expanded into 230cm recliner, with small wheels at the bottom to easy to pull and move the seat



Side V-type clamp

V-type clamp, self-locking system, 4pcs small wheels leave the ground after being locked



Castor

Equipped with 4pcs small casters for easy extension of the chair



Headrest

The design of headrest is appropriate to protect cervical vertebra and make your sleep comfortable.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.