

MK-A09

Accompany chair



*Better care
Better comfort!*

MK-A09 ACCOMPANY CHAIR

VIP concept, comfort are the main features for the MK-A09 attendant chair, it serves a chair as well as a sofa bed after stretching out.

Main Features

- Accompany chair can be open as a simple hospital bed, it can be folded into a chair.
- The structure is made of coated steel.
- Covered by PVC artificial leather, inside is high density foam.
- wooden armrest offer deluxe and gentle touch sense.
- Two pieces noiseless wheels.
- Leather color is optional.

Technical parameters

Model	MK-A09
Folding Size	L75*W80*H92cm
Unfolding Size	L185*W80*H44cm
Armrest Height	65cm
Seat width	58cm
Seat deep	48cm

Main details




Seat extension

Chair can be expanded into 185cm recliner, with small wheels at the bottom to easy to pull and move the chair



Chair position

When folded, it is a chair with tilted backrest and saves space.



Drawstring

There is a drawstring in the middle of the chair to facilitate the chair to be retracted to the seat position.



Armrest

Wide wooden armrests are more ergonomically designed and more comfortable.



Pillow

The chair has a small pillow at the head when unfolded.



Castor

With 2pcs 50 diameter casters for easy movement when unfolded.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US