

# MK-B10

## Swaddling Table



*Better care  
Better comfort!*

# MK-B10

# SWADDLING TABLE

## DESCRIPTION

The Swaddling Table MK-B10 is a specialized hospital baby swaddling station. It is designed to efficiently swaddle and organize newborns and provides a safe and comfortable space for newborns while ensuring seamless handling and care. Combining efficient swaddling capabilities with features like drawers, cabinets, and a baby dressing function can greatly enhance the care provided to newborns and make the workflow of hospital staff more streamlined.

## FEATURES

- Powder spray structure
- Removable Leather Surface allows for easy cleaning and disinfecting
- With 2 top drawers
- With 1 lockable cabinet
- Ergonomic Design provides a comfortable and practical workspace
- Guardrail with measuring tape for easy measurement
- Extension working surface increase the working area.

## STANDARD SPECIFICATION

Material	Powder spray
Drawer	2 pcs
Cabinet	1 pc
Overall size	1200*605*1050mm
Load capacity	35 KG
Mattress thickness	5cm



**Storage drawer**

The Storage Drawer provides a convenient and organized space to store all the essential items required for swaddling. With ample room, it ensures easy access to swaddle cloths, blankets, and other necessary accessories, making the swaddling process efficient and seamless.



**Storage cabinet**

The Storage Cabinet is designed to securely store the most crucial swaddling essentials. Equipped with a lock, it ensures the safety and confidentiality of these items.



**Leather mattress**

The leather surface allows for easy wiping and sanitizing, ensuring a clean and safe environment for newborns, and the mattress is removable for easy cleaning.



**EXtension working surface**

There is a extension working surface on the right side to increase the working area.



**Guardrails**

The two sides of guardrails provide a protective barrier to prevent any possible incidents of a baby falling off. Measuring tape makes it easy to measure the height of newborns.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US