

# MK-B11

Mobile Hospital Bassinet



*Better care  
Better comfort!*

# MK-B11

# MOBILE HOSPITAL BASSINET

## DESCRIPTION

The MK-B11 hospital bassinet with heavy duty frame and practical Trendelenburg position, as well as the ergonomic push handle, can bring a smoother working experience for caregivers. The full-length transparent guardrail and special double-layer guardrail design shorten the distance between mother and baby while ensuring the safety of the baby.

## FEATURES

- Powder coated steel structure
- MS perforated lying surface
- Full length of transparent guardrails around the bed
- Double-layer guardrail on the left and right sides, to ensure safety
- 10 cm double side castors with diagonal braking system
- With storage basket

## STANDARD SPECIFICATION

|                               |                          |
|-------------------------------|--------------------------|
| Overall size                  | L 1030 * W 610* H 940 mm |
| Trendelenburg                 | 0~15°±2°                 |
| Lying surface                 | L 890 * W 550 mm         |
| Height of the big guardrail   | 355 mm                   |
| Height of the small guardrail | 7 cm                     |
| Castor                        | Dia 100 mm               |



**Trendelenburg position**

Trendelenburg position can prevent baby from choking after spitting up milk.



**Gas spring lever**

Trendelenburg position is adjusted by gas spring lever.



**Transparent Guardrails**

The surrounding transparent guardrails can be lowered to facilitate nursing staff to enter and clean the bed surface.



**Double-Layer Guardrails**

The double-layer guardrail on the left and right sides facilitates caregivers to enter the crib, while the second layer of smaller, removable guardrails ensures the safety of the baby.



**Platform**

MS perforated platform is breathable and easy to clean.



**Push handle**

There are ergonomically designed push handles around the bed.



**Storage space**

Storage basket under the bed.



**Castors**

Four units of 10cm double side castors with diagonal locking system.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US