

# MK-IS13

IV stand



*Better care  
Better comfort!*

# MK-IS13

/ IV STAND



The heavy-duty base with casters and adjustable height makes the MK-IS13 suitable for various occasions in the hospital.

---

## Main Features

- Strong stainless steel base
  - Height adjustable
  - Stainless steel tubing
  - 4 foldable stainless steel hooks
  - Easy-rolling castors provide quiet effortless movement
- 

## Technical parameters

Model	MK-IS13
Height adjustment	1195~2100mm
Base size	Dia 500mm
Loading capacity	10kg
Castor	Dia 50mm
Hook	4 hooks

## Main details



### Hook

Four foldable and removable strong stainless hooks



### Knob

Height adjusted by knob



### Base

Stainless steel base with five legs , more stable and not easy to tip over.



### Castors

5 units castors with two brakes , easy to move.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US