

MK-N03

Mobile Medical Privacy Screens



*Better care
Better comfort!*



Mobile Medical Privacy Screens

Fold ▪ Easy to clean ▪ Easy to move ▪ Easy to store

MK-N03 is opening a new revolution for hospital privacy management specially in aera like patient ward, emergency room, consulting room etc..

Not only providing the privacy but also providing an enjoyable enviroment is a creative feature for this new

screen, different colorful paintings and colors are optional with customization.

Hygiene management is a superior advantage for this new model over traditional curtain due to special materials adopted which reduces tremendous contamination.

Features

- The screen can be mounted on a trolley to provide mobile privacy when necessary.
- This trolley is intended for Folding Screens in various lengths: from 5 to 11 sections.
- The sturdy frame is made of profiled aluminium tubing
- Height: 1.45 m, 1.65 m or 1.85 m.
- The patented design takes up a minimal amount of floor space, while providing maximum stability.
- Color option: White,Pink,Blue & Orange etc.

Product Details



Flexible Latch

Adjustable carrying handles can be adjusted in length and angle for maximum convenience when going up and down stairs.



Castor

The basement with four Double-sided silent castors , each castor with the brake.



Handle

The streamlined push handle is beautiful and at the same time improves the comfort when pushing and pulling.



Material

The material of the screen is aluminum profile and antibacterial plate, which is light in weight and convenient for scrubbing and disinfection.



Basement

Large "X" structure, heavy weight can make the whole screen more stable when pushed.

Optional colors



Optional patterns



Product Size

	Mobile	Wall Mounted
4P-120cm		
6P-180cm		
8P-240cm		
10P-300cm		

Technical Parameter

Mobile

Model No. Panels Only	Length	Height
(MK) N03-104(+1)	1.2M	1.45M
(MK) N03-204(+1)	1.2M	1.65M
(MK) N03-304(+1)	1.2M	1.85M
(MK) N03-106(+1)	1.8M	1.45M
(MK) N03-206(+1)	1.8M	1.65M
(MK) N03-306(+1)	1.8M	1.85M
(MK) N03-108(+1)	2.4M	1.45M
(MK) N03-208(+1)	2.4M	1.65M
(MK) N03-308(+1)	2.4M	1.85M
(MK) N03-110(+1)	3M	1.45M
(MK) N03-210(+1)	3M	1.65M
(MK) N03-310(+1)	3M	1.85M


Wall Mounted

Model No. Panels Only	Length	Height
(MK) N03-104	1.2M	1.45M
(MK) N03-204	1.2M	1.65M
(MK) N03-304	1.2M	1.85M
(MK) N03-106	1.8M	1.45M
(MK) N03-206	1.8M	1.65M
(MK) N03-306	1.8M	1.85M
(MK) N03-108	2.4M	1.45M
(MK) N03-208	2.4M	1.65M
(MK) N03-308	2.4M	1.85M
(MK) N03-110	3M	1.45M
(MK) N03-210	3M	1.65M
(MK) N03-310	3M	1.85M


Mobile



Wall mounted



Detachable design



info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.