

MK-N04

Screen



*Better care
Better comfort!*

MK-N04

SCREEN

DESCRIPTION

Stainless steel hospital ward screens are simple to erect and relocate, offering seclusion and isolation for a variety of settings, including medical and health units, hospitals, and exam rooms.

FEATURES

- Size: 1850mm*1730mm(L*H) (3 pieces)
- Stainless steel frame
- The screen has dust-proof and antibacterial function.
- Saves space when folded.
- The cloth PVC are easy to clean
- The screen can easy movement and easy installation and disassembly.

STANDARD SPECIFICATION

Model	MK-N04
Material	Stainless steel frame with waterproof PVC
Caster	50mm swivel caster
Size	1850mm*1730mm(L*H)
Panel No#	Optional(2-foldable, 3-foldable, 4-foldable...)



Material

Thickened waterproof cloth, antioxidant, and easy to clean.



Stable Bracket

Precision welding, thick and stable, not easy to dump.



Caster

Silent Casters Flexibility, easy to move.



Easy to fold when not in use, easy to store, save space

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US