

MK-P03

Emergency Cart



*Better care
Better comfort!*

MK-P03

EMERGENCY CART

DESCRIPTION

MK-P03 Hospital Crash Cart with a full range of accessories works as an ideal platform for doctor's emergency treatment. It is durable and easy for manipulation. With fast and safe access to the necessary materials and medicines.

FEATURES

- Structure: Four steel with plastic cover columns, and in ABS body
- ABS plastic top board with arc handle and protruding edge
- Side rails: 3 pcs, Stainless steel, around the top board for protection
- Four 100mm castors and two with brakes
- With five drawers: two small size, two middle size and one big size, each inner with partitions which can be organized freely, besides, each drawer with label card holder
- With Sliding side shelf, equipment holding shelf, storage box, I.V. transfusion pole, sharp container, dust basket, Power outlet and hooks, CPR board, oxygen tank holder

STANDARD SPECIFICATION

Overall size (L*W*H)	750*480*940mm/850*520*950mm
Structure	Mainly made of ABS plastic and steel-plastic
Castors	Four 100mm castors and two with brakes



Countertop

ABS engineering plastic one-step molding countertop, with handles and stainless steel guardrail.




Extended working surface

Telescopic shelf which can increase the working surface



Locking system

Central locking system with key or disposable seal



Drawer

Each drawer can be divided by different dividers, according to your needs.



Storage boxes

Big storage space, convenient to use



Garbage can

Garbage cans with different colors facilitates garbage sorting



Basket with sharp boxes

Storage for sharp needles



Power outlet

Power outlet for equipment using ,normal length 3m



CPR board

CPR board for patient rescue.



IV pole

IV pole with 4 hooks, height adjustable.



Defibrillator shelf

Angle and height adjustable, for equipment storage.



Castor

Four diameter 100mm noiseless castors with 2 brakes

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US