

# MK-P08

Anesthesia Cart



*Better care  
Better comfort!*

# MK-P08

# ANESTHESIA CART

## DESCRIPTION

MK-P08 is a trolley designed for improving efficiency and convenience to meet the general needs of the anesthesia process. A flexible divider system with removable drawer trays, and a host of organized accessories neatly arrange the medications and all related supply items.

## FEATURES

- Structure: double aluminum alloy side columns with abs body
- Expansive work area with slide-out writing surface
- Integrated and flexible drawer divider system
- Side rails: 3 pcs, abs material, around the top board for protection
- Four 100mm castors and two with brakes
- With five drawers: two small size, two middle size and one big size, each inner with partitions which can be organized freely, besides, each drawer with label card holder
- With central locker for all drawers

## STANDARD SPECIFICATION

Model	MK-P08
Overall size	760*475*950mm
Packing size	800*550*1050mm
N.W	47kg
G.W	50kg

## DRAWER SPECIFICATION

	Height	Depth	Width
Small	6.8cm	35.5cm	43cm
Middle	11cm	35.5cm	43cm
Large	22cm	35.5cm	43cm



**Drawer Divider**

Each drawer can be divided by different dividers, according to your needs.



**Tilting Box**

The tilting box can effectively classify and place different anesthesia supplies.



**Extended Work Surface**

The pull-out extension surface is convenient for medical staff to write on it.



**Sharps Box**

Sharps box safely and effectively contain and dispose of contaminated needles, syringes, and other sharps.



**Waste Bin**

Standard with flip-top waste bin, optional foot pedal waste bin.



**Push Handle**

Swallow type handle scientific and labor-saving.



**Castor**

Luxury silent 4 inch casters, 2 with brake.



**Central lock**

One key controls the locking of multiple drawers, making it convenient for professionals to manage and use them.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US