

MK-P09

Anesthesia Cart



*Better care
Better comfort!*

MK-P09

ANESTHESIA CART

DESCRIPTION

MK-P09 is designed with a synergy of function and technology to meet the unique needs of your anesthesia department. Flexible partition systems, removable drawer trays and double-door storage cabinets provide ample space for various medical equipment and supplies.

FEATURES

- Structure: Four columns made of plastic-covered steel with ABS body
- ABS plastic top board with arc handle and protruding edge
- Side rails: 3 pcs, Stainless steel, around the top board for protection
- With transparent multi-bin containers at top of trolley
- Four 100mm castors and two with brakes
- With two drawers and two cabinets: each inner of drawer with partitions which can be organized freely, each inner of cabinet with separations, besides, each drawer with label card holder

STANDARD SPECIFICATION

Model	MK-P09
Overall size	625*475*950mm
Packing size	650*550*1050mm
N.W	47kg
G.W	50kg

DRAWER SPECIFICATION

	Height	Depth	Width
Small	6.8cm	35.5cm	43cm
Middle	11cm	35.5cm	43cm

Main details



Drawer divider

Each drawer can be divided by different dividers, according to your needs.



Tilting box

The tilting box can effectively classify and place different anesthesia supplies.



Extended work surface

The pull-out extension surface is convenient for medical staff to write on it.



Waste bin

Standard with flip-top waste bin, optional foot pedal waste bin.



Basket

Reasonable use of space to place medical supplies.



Castor

Luxury silent 4 inch casters, 2 with brake.



Central lock

One key controls the locking of multiple drawers, making it convenient for professionals to manage and use them.



Double doors cabinet

A partition inside the door can be used to store small instruments.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US