

# MK-P17

Emergency Cart



*Better care  
Better comfort!*

# MK-P17

# EMERGENCY CART

## DESCRIPTION

MK-P17 crash carts is designed to deliver the ultimate in convenience and safety. This cart offers stable, smooth mobility and enhanced efficiency with multi accessories.

## FEATURES

- Made of aluminum alloy frame and abs plastic structure, light and firm
- Lightweight design and profile simplify movement and steering
- All-in-one table top with raised sides to eliminate items rolling off
- 5 abs drawers slide smooth and silent with centralized lock and plug-in labels to avoid dust
- Whisper-quiet rolling convenience
- Inner with partitions can be well organized

## STANDARD SPECIFICATION

Item	Emergency crash cart
Model	MK-p17
Overall Size (L*W*H)	850*520*990mm
Structure	Mainly made of abs plastic and steel-plastic
Castors	Four 100mm castors and two with brakes



**Countertop**

Abs handrails on the top without leakage, the sunken design is to avoid stuff falling.




**Extended Working Surface**

Extended working area with slide-out work surface.



**Locking System**

Securing all drawers right with central key locking system.



**Drawer**

5 drawers with dividers inner and card slot outside, easy for storage and identification of medicines.



**Storage Boxes**

Big storage space, convenient to use.



**Trash Can**

Trash cans with different colors facilitates trash sorting.



**CPR Board**

CPR Board provides strong, firm and stable surface.



**Basket With Sharp Boxes**

Storage for sharp needles.



**Power Outlet**

Power outlet for equipment using, normal length 3m.



**Oxygen Holder**

With belts for fixing.



**Defibrillator Shelf**

Angle and height adjustable, for equipment storage.



**Castor**

Four diameter 100mm noiseless castors with 2 brakes.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US