

MK-P28

Emergency Trolley



*Better care
Better comfort!*

MK-P28 / EMERGENCY TROLLEY

As an emergency cart, MK-P28 has all the features and flexible accessories combination. Mounting on smooth and silent casters with ergonomic push handle make it a good helper for doctors during ward rounds and emergency treatment.

Main Features

- Made of aluminum alloy frame and ABS plastic structure, light and firm
- 5 drawers with dividers inner, can be freely combined
- Detachable storage box, garbage can, basket, convenient to use
- With IV pole, Defibrillator shelf, CPR board and Oxygen tank holder
- Center-lock system with key or disposable sealing
- Four dia 100mm noiseless castors with 2 brakes
- Countertop with transparent soft glass
- Ergonomically designed handle for easy movement



Technical parameters

Model	MK-28
Overall Size	775*480*980mm
Countertop size	45*51.5cm
CPR board size	41*53cm
Extended working surface size	32*37cm
Structure	aluminum alloy frame and ABS plastic structure
Castor	Dia 100mm

Main details



Countertop

ABS engineering plastic one-step molding countertop, with guardrail and transparent soft glass.



Extended working surface

Telescopic shelf which can increase the working surface.



Locking system

Central locking system with key or disposable seal.



Drawer

5 drawers with dividers inner, easy for storage and identification of medicines.



Storage boxes

Big storage space, convenient to use.



Garbage can

Garbage cans with different colors facilitates garbage sorting.



Basket with sharp boxes

Storage for sharp needles.



Power outlet

Power outlet for equipment using, normally length 3m.



CPR board

CPR board for patient rescue.

Main details



Oxygen holder

With belts for fixing.



IV pole

IV pole with 4 hooks,height adjustable.



Defibrillator shelf

Angle and height adjustable,for equipment storage.



Castor

Four diameter 100mm noiseless castors with 2 brakes.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US