

YA-AC02

Ambulance Stretcher



*Better care
Better comfort!*

YA-AC02

Ambulance Stretcher

DESCRIPTION

The YA-AS12 ambulance stretcher is with multi functions for various environment operation. Optimized structure provides the stretcher for great stability and excellent light weight for easy operation. Auto-kick foldable design makes the stretcher operational for one medical staff from group upto ambulance automatically.

FEATURES

- Lightweight, durable aluminum construction.
- Easy-to-use release handle design
- The surface is made of foamed cushion, vinyl protective cover that is durable and easy to clean
- Swing-down siderails enable convenient patient transfer from bed to cot
- The upper platform can be easily separated
- Equip I.V. pole
- Two stationary and two 360°swivel wheels

STANDARD SPECIFICATION

Material	High-strength aluminum alloy
Size(Low position)	1950*570*280mm
Size(High position)	1950*570*910mm
Backrest adjustable	0-75°
Footrest adjustable	0-15°
Load capacity	200kg



BACKREST ADJUST BY GAS SPRING

The backrest is controlled by gas spring, and the back can be adjusted to the desired angle by just pressing lightly.



LEGREST ADJUST

The legrest is adjusted by the rack, and there are 2 gears to adjust the calf plate.



SIDE RAIL

The foldable side rail can prevent the patient from falling off the stretcher and causing secondary injury.



SEPARATE HANDLE

There is a red handle at the end of the stretcher, pulling the red handle can separate the upper platform from the stretcher.



LEG FOLDED HANDLE

Pressing the red handle at the same time can make the stretcher into reverse trendelenburg position; when pressing the blue handle, the legs can be quickly folded and into the ambulance.



MATTRESS

The concave mattress provides support to the patient while minimizing the shaking of the patient and avoiding secondary injuries.



CASTOR

150mm diameter casters with single brake.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US