

ABS BEDSIDE CABINET


DRAWER

The pull-out drawer is sturdy, rust-proof and easy to clean, it can hold small items



TOWEL HOOK

Hidden type towel hook and rack, firm, anti-rust, easy to clean. It can be put away when not in use



STORAGE

The cabinet is equipped with grooves for holding hot water bottles and adjustable shelves to increase storage space.



HIDDEN DINING TABLE

The hidden dining table has grooves for placing water cups and chopsticks to increase the space for patients. ABS plastic is easy to clean.



CASTOR

The bedside table has wheels for easy movement and two brakes. It can be fixed after moving to the designated position.



ABS Bedside cabinet

Model: YA-B06

The YA-B06 bedside cabinet is made of 100% new ABS plastic, equipped with a retractable dining table. It can be moved to anywhere with optional wheels. Convenience and high space utilization rates solve the problem of limited hospital space.

- Whole ABS forming.
- Cabinet top, drawer, dining table, layer board, door and base are made of ABS.
- Towel hook and rack on the 2 sides of the cabinet, that can be put away.
- Dining board is between the drawer and tabletop.
- Thermos, washbasin and other commodities can be put in it.
- With four castors, easy to move.

TECHNICAL DATA

ITEM	PARAMETERS
Size	L475*W475*H800mm
Accessories	1 drawer, 1 cabinet, 1 sliding dining table
Caster	Dia.50mm



Zhangjiagang Medi Medical Equipment Co., Ltd
 Tel: +86 15250358182
 E-mail: info@medik-medical.com

