

YA-B15

Bedside cabinet



*Better care
Better comfort!*

YA-B15 / BEDSIDE CABINET

The YA-B15 is constructed with metal frame and covered with stainless steel top for easy cleaning and heavy duty use. And spacious storage with all the features make it as an ideal option for the hospital.

Main Features

- Frame: Metal material
- Cabinet top is made of stainless steel and others are made of epoxy coating
- Thermos , washbasin and other commodity can be put in it

Technical parameters

Model	YA-B15
Overall size	500*500*800mm
Material	Stainless steel and metal epoxy coating
Configuration	1 drawer and 1 cabinet
Load capacity	50kg



Main details



Countertop

Stainless steel table top, strong and not easy to rust.



Drawer

The drawers with rails and key locks can store some small item.



Storage space

There is a shelf to increase the storage space in the cabinet.



Castor

Dia 50mm casters make it easy to move the bedside cabinet, and two of them have brakes to fix bedside cabinet.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US