

YA-D3-2

Electric Hospital Bed
With Manual Backup



*Better care
Better comfort!*

YA-D3-2

ELECTRIC HOSPITAL BED WITH MANUAL BACKUP

DESCRIPTION

YA-D3-2 is an electric hospital bed with manual backup, the special system provides double security in case of power failure or motor damage. So it is an ideal option specially for those area with unstable power supply.

FEATURES

- Ms perforated platform with ventilation hole, molded
- Full length split side rails with safety gap
- Manual crank as backup for back,foot section and height adjustment
- Auto-counter
- Detachable head/foot board with easy lock
- Central brake system with wheel diameter of 125mm
- Urine bag hooks on both sides
- Bumpers on 4 corners for protection
- Adjustable iv pole


STANDARD SPECIFICATION

Overall size	2230*1100*520-770mm
Lying size	2045*900mm
Back section up/down	0-75°
Foot section up/down	0-40°
Loading capacity	205kg



Hand remote

Hand remote for back,foot section and height adjustment as well as auto-counter function.



Manual crank

Removable crank can be inserted into any hole to control the back,foot section and height.



Urine bag hook

Urine bag hooks on both sides of the bed for easy use in any position.



Side rail

Full length split side rail, tuck away type ,with safety gap to protect patient.



IV pole

Adjustable IV pole with 4 hooks.



Bed Frame

MS perforated platform with ventilation hole,antibacterial, corrosion resistant.



Removable head and footboard

Removable head and footboard with locking mechanism.



Castors

4 luxurious silent double-sided wheels with central brake system.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US