

# YA-D5-10

Electric Hospital Bed



*Better care  
Better comfort!*

# YA-D5-10

# Electric Hospital Bed

## DESCRIPTION

The five-function electric hospital bed with soft connection, and the adjustment part has a soft connection to improve the comfort of the patient. Central lock, the detachable headboard and footboard of the bed greatly provide convenience for medical personnel.

## FEATURES

- Four split side rail, drop down under the bed platform
- Hand remote control
- Detachable head/footboard with easy lock
- With articulated ABS soft joint between different platform sections
- Central brake system with a wheel diameter of 125mm
- Auto-counter to provide the patient from sliding onto foot end
- Manual CPR on both sides(optional)

## STANDARD SPECIFICATION

Overall size	2160*1050 *450-700mm
Lying size	1910*900 mm
Backrest adjustable	0-75°
Footrest adjustable	0-40°
Height adjustable	450-700mm
Trendeleburg	0-12°
Revere- Trendelenburg	0-12°
Loading capacity	205kg



**REMOTE CONTROL**

The integrated remote control with clear icons allows nurses to quickly and clearly understand the sections they need to control.



**MOTORS**

CE mark motors can be replaced according to requirements. Strong stability, great support without any noise. Durable and high safety motor system with compliance to IEC60601 standard, low noise.



**DETACHABLE HEADBOARD**

The lock structure at the end of the bed is convenient for quick removal in an emergency, allowing medical personnel to be closer to patient's head carry out rescue.



**SOFT LINK**

Soft link makes platform more smooth to let the patient feel more comfortable when adjust the position.



**PP SAFETY RAILING**

Tuck-away PP safety railing with damper for slow dropping, special patented design to avoid accidental falling from patient's on bed operation.



**CENTRAL LOCKING SYSTEM**

Central locking system with dual sided wheels is with excellent stability and convenient for medical staff to manipulate.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 199 6269 8550**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US