

YA-D5-7

ICU bed



*Better care
Better comfort!*

YA-D5-7

ICU bed

DESCRIPTION


The YA-D5-9 ICU bed is compactly designed with elegant appearance, easy to operate. Wide range height adjustment facilitate medical operation.

FEATURES

- Cold rolled steel with powder coating treatment
- Remote control
- Detachable ABS head/foot board
- Split Guard Rails: Guard rail PP- 4pcs provides protection to patient against fall.
- Manual and electric CPR
- Degree indicator on side rails display the angle of back rest and Trendelenburg / Reverse Trendelenburg.
- Noiseless Luxury castors- Dia 125mm with Central braking system.
- Provision for I.V. Rod on all sides of the bed.

STANDARD SPECIFICATION

Overall size	2045*1060*500mm
Height adjustable	500-780mm
Backrest	0-75°
Footrest	0-45°
Trendelenburg	12°
Reverse- Trendelenburg	12°
Caster	Dia125mm
Loading capacity	250kg



Remote nurse control, electric CPR included.



Degree indicator, easy for position indication.



Detachable ABS head/foot board for emergency application.



4 pcs dia 125mm castors, central locking.



Manual CPR in each side of back board.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US