

YA-DS-R04

Electric Lifting Recliner Sofa



*Better care
Better comfort!*

YA-DS-R04

ELECTRIC LIFTING RECLINER SOFA

DESCRIPTION

Electrically lifting function to help the people from a sitting to standing position is the great feature of this recliner sofa chair via a single button, The comfortable and luxurious seat cushion is equipped with a fully electric remote control operating system and a central control caster system, making the YA-DS-R04 a luxurious Lifting Recliner Sofa suitable for long-term patient rehabilitation and care.

FEATURES

- Luxurious upholstery for a comfortable experience
- Fully electric operation, convenient for caregivers and patients
- Casters and push handles for easy mobility
- Five-position: lifted, seated, TV recliner, full recliner, zero gravity
- Lift mechanism can assist patients in standing
- Hand remote with a charging port for phone or iPad

STANDARD SPECIFICATION

Overall size(L*W*H)	1025*720*1075mm
Seat size(W*D*H)	505*535*515mm
Back height	780mm
Full recliner length	1660mm
Back section adjustment	175°
Foot section adjustment	90°
Lifted angle	45°
Full recliner length	1660mm



Upholstery

Luxurious cushions bring users more comfort.



Hand remote

Easy-to-read buttons and a charging port for charging your phone or iPad.



Push handle

The back of the sofa is equipped with a push handle for easy pushing.



Brake

Lock the casters by pressing and releasing the crossbar.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.