

YA-ET-D01

Electric Examination Treatment Table



*Better care
Better comfort!*

YA-ET-D01 / EXAMINATION TREATMENT TABLE

Comfort and functionality have always been factors that must be considered to choose a proper treatment table. As an economical electric treatment table, YA-ET-D01 has the necessary height and back adjustment functions, is equipped with luxurious and comfortable Upholstery, heavy duty structure and a stable motor system, and is suitable for various treatments, examinations, and salons.

Main Features

- Features a 2-section top design
- Electric height and back section adjustment with hand remote
- Face hole allows for easy prone position
- Ultra-strong, uni-frame base and top construction
- Seamless, rounded corner, upholstered top
- System of retractable castors that are activated through a pedal

Technical parameters

| | |
|------------------|---------------------|
| Model | YA-ET-D01 |
| Overall size | 184*68*48~85 cm |
| Back section | 0-65°±5° |
| Castors | Dia 75mm castors *2 |
| Loading capacity | 180 kg |

Function



• Height adjustment



• Back section adjustment

Main details



Upholstery

Luxurious and seamless upholstery to improve the user's comfort and for easy clean.



Face Hole

Face hole for prone treatment.



Hand remote

Height and back section adjusted by hand remote.



Base

2 front retractable wheels and 4 feet anchored to the floor, easy to move.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US