

# YA-H2-1

Manual Bed



*Better care  
Better comfort!*

# YA-H2-1

# MANUAL BED

## DESCRIPTION

The YA-H2-1 Full Fowler bed with an ergonomic and compact design, it can meet the general demand of hospitalization, bed frame with electrostatic epoxy coating can prevent it from corrosion. mesh platform provides good ventilation for patients and is more suitable for long-term bedridden patients.

## FEATURES

1. Bed Frame electrophoresis and powder coating
2. Mesh bed platform, easy to ventilate
3. Stainless steel with MDF Headboard and footboard, removable
4. The back section adjusted by ratchet
5. Height adjusted by crank
6. Castors 125mm diameter, 1 guide and 2 with brake

## STANDARD SPECIFICATION

Overall Size	2100*950*500-800mm
Lying surface	1950*900mm
Back section adjustable	0° to 60°
Height adjustable	500-800mm
Safe working load	200kg



**Mesh platform**

Mesh platform provides great ventilation features to the patient especially when a long time in bed is needed.



**Back section**

The back section is adjusted by ratchet and has multiple gears.



**Crank**

The metal foldable crank can prevent the handle from protruding outside to avoid collision.



**Gas spring**

The crank has been enhanced with a gas spring to make the bed easier and smoother when raising or lowering its height.



**Castor**

Dia 125mm castors, 2pcs with brake and 1pc with guide function.



**Head/foot board**

stainless steel head/foot board with laminated panel gives an attractive appearance and easy to push the bed.



**IV pole hole**

There are 4 infusion holes around the bed, which is convenient for placing the IV pole.



**Bumper castor**

Bumper castors on the four corners of the bed, which can reduce the impact of accidentally hitting the wall during the movement.



**IV pole**

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US