

YA-M1-5

Manual hospital bed



*Better care
Better comfort!*

YA-M1-5

MANUAL HOSPITAL BED

DESCRIPTION

Electrostatic epoxy coating finishing, durable foldable crank and safe guardrail make the YA-M1-5 as a basic general ward bed with low budget.

FEATURES

- Bed Frame electrophoresis and powder coating.
- One-off forming of carbon steel bed plate with large airhole design, well ventilated and great strength.
- ABS head and foot boards formed by blow moulding.
- Back section adjusted by crank.
- Aluminium alloy foldable side rails are durable and easy to operate.

STANDARD SPECIFICATION

Overall Size	2000*900*530mm
Back section adjustable	0° to 75°
Fixed height	530mm
Safe working load	180kg



Slatted platform

Die punching bed board features a robot-welded surface that has been polished to perfection, leaving no welding scars.



crank

The metal foldable crank can prevent the handle from protruding outside to avoid collision.



Leg

Four legs anchored to the floor with anti-slip foot cover



Head/foot board

ABS head and foot boards formed by blow moulding, the head and footboards are effortless to disassemble and clean.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US