

medik

Tel : +86 15250358182

E-mail : info@medik-medical.com




YA-M3-3

MANUAL HOSPITAL BED




MAIN FEATURES

- * Backrest, Footrest and Trendelenburg movements by 3 cranks
- * Epoxy coated mild steel bed frame
- * 4 section bed board
- * ABS side rails with angle indicator
- * ABS head and foot board with holes for moving the bed easily
- * Dia125mm noiseless castors ,all with individual brakes
- * Bumper castors on 4 corners for protection



BACK SECTION



FOOT SECTION



HEIGHT ADJUSTMENT

PRODUCT DETAILS



FOLDABLE CRANK HANDLE



CASTORS



HEADBOARD / FOOTBOARD



BUMPER COVER



ADJUSTING SYSTEM



STEEL | PAINTING

TECHNICAL PARAMETER

ITEM	PARAMETER
Model	YA-M3-3
Overall size	2160*960*500mm
Back section adjustable	0-75°
Foot section adjustable	0-40°
Trendelenburg	0-12°
Castor	Dia125mm
Safe working load	175kg