

YA-M3-4

Manual Hospital Bed



*Better care
Better comfort!*

YA-M3-4

Manual Hospital Bed

DESCRIPTION

The YA-M3-4 is a manual hospital bed with ergonomic design and has elegant appearance, 4units PP side railing with angle display provides the safety for patients and easy for bed exit, it is widely accepted by hospital ward area with limited budget.

FEATURES

- Bed Frame electrophoresis and powder coating
- MS perforated platform with ventilation hole, molded
- Four split side rail, drop down under the bed platform
- Head and Footboard Detachable head/footboard with easy lock
- 3 pieces foldable hand cranks
- 6 urine bag hooks on both sides
- With angle indicator embedded in the railing for back and foot

STANDARD SPECIFICATION

Model	YA-M3-4
Overall size	2240*1085*450-700mm
Recommended mattress size	2015*875mm
Height adjustable	450-700mm
Backrest adjustable	0° to 70°
Footrest adjustable	0° to 40°
Loading capacity	200kg



ANGLE DISPLAY

The side rail with 3 gears to adjust the height, anti-pinch hands to protect children's safety . When adjusting, It needs to be pulled two handles at the same time to adjust the side rail which prevents children from misoperation.



IV HOLE

There are infusion holes on the four corners of the bed, and two more at the waist position, which can provide more solutions.



URINE BAG HOOK

There are 6 urine bag hooks on both sides of the bed for easy use in any position.



SIDE RAIL LOCK

Tuck-away side railing meets IEC60601-2-52 safety standard, special protection mechanism can prevent patient to drop down the railing from inside bed.



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 199 6269 8550

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US