

YA-M5-1

Five Function Manual Hospital Bed



*Better care
Better comfort!*

YA-M5-1

Five function Manual hospital bed

DESCRIPTION

YA-M5-1 is a five-crank manual hospital bed with an ergonomic design and modern design. It has all the functions required for medical treatment. 4 pieces of full-length ABS side rails with angle display provide safety and convenience for patients to exit the bed. The removable head and foot panels and central control casters make it more efficient and flexible during the process of care.

FEATURES

- MS perforated platform with ventilation hole, molded
- Four split side rails, drop down under the bed platform
- 4 pieces of foldable hand cranks
- Detachable head/foot board with easy lock
- Central brake system with wheel diameter of 125mm
- Urine bag hooks on both sides
- With angle indicator embedded in the railing
- Bumper castors on 4 corners for protection
- Retractable 5th wheel for direction(optional)

STANDARD SPECIFICATION

Overall size	2490*1100 *520-770mm
Lying size	2300*900 mm
Back section up/down	0-75°
Foot section up/down	0-40°
Height adjustment	520-770mm
Trendelenburg	0-12°
Revere- Trendelenburg	0-12°
Loading capacity	195kg



Manual crank

Four pieces of foldable manual crank at foot side of the bed.



Angle display

Angle display on the side rail, can indicate the adjustment angle of the bed.



Urine bag hook

Urine bag hooks on both sides of the bed for easy use in any position.



Side rail

Four split side rails, tuck-away type, special protection mechanism can prevent patient to drop down the railing from inside bed.



Bumper castors

Bumper castors on 4 corners for protection.



Bed Frame

MS perforated platform with ventilation hole, antibacterial, corrosion resistant.



Castors

4 luxurious silent double-sided wheels with central brake system.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US