


YA-M5-6

Manual hospital bed



*Better care
Better comfort!*

YA-M5-6

MANUAL HOSPITAL BED

DESCRIPTION

The YA-M5-6 is a 5 functions manual hospital bed with ergonomic design elegant appearance. Patient safety and functionality at low cost are the priority concerns for this bed. It provides a complete range of bed postures with low budget.

FEATURES

- Bed Frame electrophoresis and powder coating
- MS perforated platform with ventilation hole, molded
- Four split side rail, drop down under the bed platform
- Head and Footboard Detachable head/footboard with easy lock
- 4 pieces foldable hand cranks
- urine bag hooks on both sides
- With angle indicator embeded in the railing
- Oxygen bottle holder and IV pole

STANDARD SPECIFICATION

Overall size	2240*1085*450-700mm
Recommended mattress size	2015*875mm
Height adjustable	450-700mm
Backrest adjustable	0° to 70°
Footrest adjustable	0° to 40°
Trendelenburg/reverse trendelenburg	0-12°
Castor	Dia 125mm
Loading capacity	200kg



Section adjustable

Back, foot section and height can be adjusted by crank separately, can meet the needs of patients in different positions.



Trendelenburg/reverse trendelenburg

Trendelenburg/reverse trendelenburg are adjusted by manual crank.



Manual crank

Four pieces foldable manual crank at foot side of the bed, and there is a medical record holder at the end of the bed for easy identification by medical staff.



Angle display

Back, foot section and height can be adjusted by crank separately, can meet the needs of patients in different positions.




IV hole

There are infusion holes on the four corners of the bed.



Urine bag hook

Urine bag hooks on both sides of the bed for easy use in any position.



Oxygen bottle holder

There is an oxygen bottle holder at the bottom of the bed for easy use.



Side rail lock

Tuck-away side railing meets IEC60601-2-52 safety standard, special protection mechanism can prevent patient to drop down the railing from inside bed.



Castors

4 luxurious silent double-sided wheels with individual brake.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US