

YA-PM1-1

Pediatric bed



*Better care
Better comfort!*

YA-PM1-1

Pediatric bed

FEATURES

- Cold rolled steel with powder coating treatment
- Back board adjustment by manual crank
- 4 pcs swivel type castors, all with locking
- Platform with ventilator hole design
- Full length side railing for safety protection

STANDARD SPECIFICATION

Overall size	1680*870*500mm
Platform size	1550*800mm
Back board	0-75°
Height (from floor to top)	1100mm
Height (from platform to top)	600mm
Castor size	125mm
Safety loading weight	100kg



Foldable manual crank for back board adjustment.



4 pcs dia 125mm castors, swivel type with locking.



Cold rolled steel platform with ventilator hole design, easy for skin breath.



Safety full length side railing will protect patient well from falling down, easy for up-down.



Head board is made of anti-bacterial HPL board, easy for clean and sterilization.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US