

YA-PM1-2

Pediatric Bed



*Better care
Better comfort!*

YA-PM1-2

PEDIATRIC BED

DESCRIPTION

YA-PM1-2 is a pediatric bed with 360°protected guardrail for great security. It is a very durable and pediatric hospital ward bed with smooth backrest adjustment.

FEATURES

- Bed frame made of stainless steel materials.
- Back board adjustment by manual crank.
- 4 pcs swivel type castors, two with locking.
- Platform made of cold-rolled steel plate.
- Full length side railing for safety protection.

STANDARD SPECIFICATION

Overall size	1600*1000*700mm
Railing height	530mm
Back adjustment	0-75°
Castor size	φ75mm
Loading capacity	80kg



One chrome plated manual crank for back board adjustment



4 pcs dia 75mm castors, two with locking



Height-adjustable IV pole with 4 foldable stainless steel hooks, easy for managing medical equipment

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US