

YA-PM3-3

Pediatric Bed



*Better care
Better comfort!*

YA-PM3-3

PEDIATRIC BED

DESCRIPTION

The YA-PM3-3 pediatric bed with trendelenburg and reverse trendelenburg function bears all the necessary features for a standard Pediatric Bed, and the high full-length siderail is designed to prevent children from falling. Comfortably for children's recovery, economically for hospitals to buy.

FEATURES

- Bed Frame electrophoresis and powder coating
- Bed platform MS perforated platform with ventilation hole, molded
- Removable MDF Headboard and footboard, cartoon print
- Full-length side rail can be adjusted up and down
- Back section adjusted by ratchet
- Trendelenburg and reverse Trendelenburg adjust by crank
- Castors 150mm diameter with central braking system

STANDARD SPECIFICATION

Overall Size	1800*960*600mm
Back section adjustable	0° to 60°
Trend&reverse trend	0°to 10°
Safe working load	150kg



Bed platform

MS perforated bed platform with ventilation hole, molded.



crank

Foldable metal cranks for trendelenburg and reverse trendelenburg adjustments.



Castor

Dia 150mm castors with central braking system.



Side rail

The full-length side rail can be adjusted up and down to prevent children from falling.



Head/foot board

stainless steel head/foot tube with cartoon painting panel gives an attractive appearance and easy to push the bed.



IV stand



Monitor Tray



Oxygen cylinder holder

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US