

YA-ST03

Hydraulic Shower Trolley



*Better care
Better comfort!*

YA-ST03

hydraulic shower trolley

DESCRIPTION

How to ensure the safety of elderly people and patients with limited mobility during bathing has always been a matter of great concern to all medical staff. Now, YA-ST03 shower trolley provides a solution, the stable double-column hydraulic adjustment system, all-round guardrail protection and easy transfer system can bring more efficient and safe bathing care to patients and medical staff.

FEATURES

- Full length stainless steel guardrail, rust-proof and high safety.
- Dual hydraulic control, stable and convenient.
- Waterproof bed surface, plugs and drain pipes for convenient bathing care.
- With trendelenburg and reverse trendelenburg function to promote drainage.
- 200mm diameter wheels, convenient for transportation.
- 5th castor with steering function.

STANDARD SPECIFICATION

Overall size	L214*W83*H73~105±3cm
Platform size	L196*W60*H57~88±3cm
Trend/ rev	-12°~ 12°±2°
Castor	Dia 200mm
Safe working load	240KG



Soft mattress

The soft mattress with waterproof cover can realize bathing function.



Plug and drain hose

Equipped with the plug and drain hose, can facilitate drainage.



Guardrail

Full-length stainless steel guardrail all around the trolley, can be raised and dropped down, rust-proof and provide safety.



Guardrail adjustment device

The guardrails around the bed can be easily raised or lowered by mechanical devices, also with locking function.



Foot pedal

The foot pedals around the bed for easy access.



Fifth castor

The fifth wheel has a guiding function.



Castors

Diameter 200mm castors design with central locking and steering facility.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US